

Outer membrane protein 26 H. pylori protein (STJ504237)

STJ504237

GENERAL INFORMATION

Product Type	Proteins
Short Description	Outer membrane protein 26 H. pylori-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Heat killed and denatured Helicobacter pylori extracts
Immunogen Region	
Specificity	
Immunogen Sequence	