

CMT protein (STJ504241)

STJ504241

GENERAL INFORMATION

Product Type	Proteins
Short Description	CMT-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23463
Gene Symbol	ICMT
Uniprot ID	ICMT_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 230-284 on human Protein-S-isoprenylcysteine O-methyltransferase.
Immunogen Region	230-284
Specificity	
Immunogen Sequence	