

IL12 p40 protein (STJ504248)

STJ504248

GENERAL INFORMATION

Product Type	Proteins
Short Description	IL12 p40-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region spanning 280-335 on rat Interleukin-12 subunit beta protein.
Immunogen Region	280-335
Specificity	
Immunogen Sequence	