

IL17RB protein (STJ504251)

STJ504251

GENERAL INFORMATION

Product Type	Proteins
Short Description	IL17RB-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55540
Gene Symbol	IL17RB
Uniprot ID	I17RB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-400 on IL17RB protein.
Immunogen Region	300-400
Specificity	
Immunogen Sequence	