

IL17C protein (STJ504252)

STJ504252

GENERAL INFORMATION

Product Type	Proteins
Short Description	IL17C-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84818
Gene Symbol	IL17RC
Uniprot ID	I17RC_HUMAN
Immunogen	Synthetic peptide corresponding to IL17RC protein common to C1, C2, C5, C8, C9, C11, C12 isoforms.
Immunogen Region	
Specificity	
Immunogen Sequence	