

IL23 P19 protein (STJ504255)

STJ504255

GENERAL INFORMATION

Product Type	Proteins
Short Description	IL23 P19-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	51561
Gene Symbol	IL23A
Uniprot ID	IL23A_HUMAN
Immunogen	Synthetic peptide corresponding to C terminal region amino acids 162-179 of Human IL23 P19. Peptide was covalently modified. Sequence: FKILRSLQAFAAVAARV
Immunogen Region	162-179
Specificity	
Immunogen Sequence	