

Integrin linked ILK protein (STJ504257)

STJ504257

GENERAL INFORMATION

Product Type	Proteins
Short Description	Integrin linked ILK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3611
Gene Symbol	ILK
Uniprot ID	ILK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-452 on human Integrin-linked protein kinase.
Immunogen Region	400-452
Specificity	
Immunogen Sequence	