

inaD protein (STJ504258)
STJ504258**GENERAL INFORMATION**

Product Type	Proteins
Short Description	inaD-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to C terminal amino acids 659-674 of Fruit fly (<i>Drosophila melanogaster</i>) inaD protein. Sequence: EVTRPKPTLRTEAPKA
Immunogen	659-674
Region	
Specificity	
Immunogen	
Sequence	