

Inhibin alpha protein (STJ504259)

STJ504259

GENERAL INFORMATION

Product Type	Proteins
Short Description	Inhibin alpha-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3623
Gene Symbol	INH A
Uniprot ID	INH A_HUMAN
Immunogen	N-epitope synthetic peptide taken within amino acid region 1-50 on Inhibin alpha gene.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	