

iNOS protein (STJ504264)

STJ504264

GENERAL INFORMATION

Product Type	Proteins
Short Description	iNOS-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4843
Gene Symbol	NOS2
Uniprot ID	NOS2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1064-1114 on human inducible nitric oxide synthase protein.
Immunogen Region	1064-1114
Specificity	
Immunogen Sequence	