

Integrin beta 3 protein (STJ504266)

STJ504266

GENERAL INFORMATION

Product Type	Proteins
Short Description	Integrin beta 3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3690
Gene Symbol	ITGB3
Uniprot ID	ITB3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human Integrin Beta 3 subunit protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	