

IRE1A protein (STJ504268)

STJ504268

GENERAL INFORMATION

Product Type	Proteins
Short Description	IRE1A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2081
Gene Symbol	ERN1
Uniprot ID	ERN1_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence within amino acid region 900-1000 on human IRE1A protein.
Immunogen Region	900-1000
Specificity	
Immunogen Sequence	