

ITK protein (STJ504271)

STJ504271

GENERAL INFORMATION

Product Type	Proteins
Short Description	ITK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3702
Gene Symbol	ITK
Uniprot ID	ITK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 550-600 on human ITK protein.
Immunogen Region	550-600
Specificity	
Immunogen Sequence	