

ITPR2 protein (STJ504273)

STJ504273

GENERAL INFORMATION

Product Type	Proteins
Short Description	ITPR2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81678
Gene Symbol	Itpr2
Uniprot ID	ITPR2_RAT
Immunogen	Synthetic peptide taken within amino acid region 100-150 on rat Inositol 1, 4, 5-triphosphate receptor type 2 protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	