

KDM4C/GASC1/JMJD2C protein (STJ504275)

STJ504275

GENERAL INFORMATION

Product Type	Proteins
Short Description	KDM4C/GASC1/JMJD2C-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23081
Gene Symbol	KDM4C
Uniprot ID	KDM4C_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-520 on human JmjC domain-containing histone demethylation protein 3C protein.
Immunogen Region	500-520
Specificity	
Immunogen Sequence	