

JMJD2D protein (STJ504276)

STJ504276

GENERAL INFORMATION

Product Type	Proteins
Short Description	JMJD2D-protein if for use in research applications.
Applications	Immuno-depletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55693
Gene Symbol	KDM4D
Uniprot ID	KDM4D_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 340-390 on human Jumonji domain containing 2D protein.
Immunogen Region	340-390
Specificity	
Immunogen Sequence	