

Ki67 protein (STJ504278)

STJ504278

GENERAL INFORMATION

Product Type	Proteins
Short Description	Ki67-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4288
Gene Symbol	MKI67
Uniprot ID	KI67_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 2950-3000 on human antigen KI-67 isoform 1.
Immunogen Region	2950-3000
Specificity	
Immunogen Sequence	