

KIF13A protein (STJ504279)

STJ504279

GENERAL INFORMATION

Product Type	Proteins
Short Description	KIF13A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	63971
Gene Symbol	KIF13A
Uniprot ID	KI13A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1700-1800 on human KIF13A isoform a protein.
Immunogen Region	1700-1800
Specificity	
Immunogen Sequence	