

KIF14 protein (STJ504280)

STJ504280

GENERAL INFORMATION

Product Type	Proteins
Short Description	KIF14-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9928
Gene Symbol	KIF14
Uniprot ID	KIF14_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 850-960 on human KIF14 protein.
Immunogen Region	850-960
Specificity	
Immunogen Sequence	