

p27 KIP1 protein (STJ504281)

STJ504281

GENERAL INFORMATION

Product Type	Proteins
Short Description	p27 KIP1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	N-epitope synthetic peptide corresponding to unique amino acid sequences on human p27 KIP1 protein.
Immunogen Region	N-Term
Specificity	
Immunogen Sequence	