

LAMP2 protein (STJ504289)

STJ504289

GENERAL INFORMATION

Product Type	Proteins
Short Description	LAMP2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:50-1:100
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3920
Gene Symbol	LAMP2
Uniprot ID	LAMP2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 260-310 on human lysosomal membrane glycoprotein-2.
Immunogen Region	260-310
Specificity	
Immunogen Sequence	