

LAT protein (STJ504292)

STJ504292

GENERAL INFORMATION

Product Type	Proteins
Short Description	LAT-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27040
Gene Symbol	LAT
Uniprot ID	LAT_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 40-90 on LAT protein.
Immunogen Region	40-90
Specificity	
Immunogen Sequence	