

GPCR LGR6 protein (STJ504294)

STJ504294

GENERAL INFORMATION

Product Type	Proteins
Short Description	GPCR LGR6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	59352
Gene Symbol	LGR6
Uniprot ID	LGR6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 850-920 on human LRG6 protein.
Immunogen Region	850-920
Specificity	
Immunogen Sequence	