

LILRB3 protein (STJ504298)

STJ504298

GENERAL INFORMATION

Product Type	Proteins
Short Description	LILRB3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	102725035/107987462/11025
Gene Symbol	LILRB3
Uniprot ID	LILRB3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 580-631 on human LILRB3 protein.
Immunogen Region	580-631
Specificity	
Immunogen Sequence	