

LOXL1 protein (STJ504302)

STJ504302

GENERAL INFORMATION

| | |
|--------------------------|----------------------------------------------------|
| Product Type | Proteins |
| Short Description | LOXL1-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|-----------------------------------------------------------------|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--------------------------------------------------------------------------------|
| Gene ID | 16949 |
| Gene Symbol | Lox1 |
| Uniprot ID | LOXL1_MOUSE |
| Immunogen | Synthetic peptide taken within amino acid region 570-607 on rat LOXL1 protein. |
| Immunogen Region | 570-607 |
| Specificity | |
| Immunogen Sequence | |