

LARP7 protein (STJ504343)

STJ504343

GENERAL INFORMATION

Product Type	Proteins
Short Description	LARP7-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51574
Gene Symbol	LARP7
Uniprot ID	LARP7_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	