

NTAL protein (STJ504345)
STJ504345**GENERAL INFORMATION**

Product Type	Proteins
Short Description	NTAL-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7462
Gene Symbol	LAT2
Uniprot ID	NTAL_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	