

Alpha Fetoprotein protein (STJ504357)

STJ504357

GENERAL INFORMATION

Product Type	Proteins
Short Description	Alpha Fetoprotein-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11657
Gene Symbol	Alb
Uniprot ID	ALBU_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 50-100 on mouse alpha-fetoprotein.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	