

MAP3K10 protein (STJ504358)

STJ504358

GENERAL INFORMATION

Product Type	Proteins
Short Description	MAP3K10-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4294
Gene Symbol	MAP3K10
Uniprot ID	M3K10_HUMAN
Immunogen	Synthetic peptide corresponding to the C-terminal taken within amino acid region 900-950 on human MAP3K10 protein.
Immunogen Region	900-950
Specificity	
Immunogen Sequence	