

MAP3K14 protein (STJ504360)

STJ504360

GENERAL INFORMATION

Product Type	Proteins
Short Description	MAP3K14-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9020
Gene Symbol	MAP3K14
Uniprot ID	M3K14_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-350 on human MAP3K14 protein.
Immunogen Region	250-350
Specificity	
Immunogen Sequence	