

MC3 Receptor protein (STJ504364)

STJ504364

GENERAL INFORMATION

Product Type	Proteins
Short Description	MC3 Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4157
Gene Symbol	MC1R
Uniprot ID	MSHR_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequences on MC3r protein.
Immunogen Region	
Specificity	
Immunogen Sequence	