

Mitofusin 2 protein (STJ504372)

STJ504372

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | Mitofusin 2-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 9927 |
| Gene Symbol | MFN2 |
| Uniprot ID | MFN2_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 543-593 on human Mitofusin 2 protein. |
| Immunogen Region | 543-593 |
| Specificity | |
| Immunogen Sequence | |