

MERTK protein (STJ504375)

STJ504375

GENERAL INFORMATION

Product Type	Proteins
Short Description	MERTK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4507
Gene Symbol	MTAP
Uniprot ID	MTAP_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-200 on MerTK protein.
Immunogen Region	100-200
Specificity	
Immunogen Sequence	