

Myeloperoxidase protein (STJ504379)

STJ504379

GENERAL INFORMATION

Product Type	Proteins
Short Description	Myeloperoxidase-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4353
Gene Symbol	MPO
Uniprot ID	PERM_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 600-650-625 human Myeloperoxidase protein
Immunogen Region	600-650-625
Specificity	
Immunogen Sequence	