

MRGX1 protein (STJ504384)

STJ504384

GENERAL INFORMATION

Product Type	Proteins
Short Description	MRGX1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	282547/404646
Gene Symbol	Mrgrpx1
Uniprot ID	MRGX1_RAT
Immunogen	Synthetic peptide taken within amino acid region 280-330 on rat MrGX1 protein.
Immunogen Region	280-330
Specificity	
Immunogen Sequence	