

MUC1 protein (STJ504390)

STJ504390

GENERAL INFORMATION

Product Type	Proteins
Short Description	MUC1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4582
Gene Symbol	MUC1
Uniprot ID	MUC1_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on MUCIN1 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	