

Myoglobin protein (STJ504393)

STJ504393

GENERAL INFORMATION

Product Type	Proteins
Short Description	Myoglobin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4151
Gene Symbol	MB
Uniprot ID	MYG_HUMAN
Immunogen	Full length affinity purified human myoglobin protein.
Immunogen Region	Full length
Specificity	
Immunogen Sequence	