

Myocilin protein (STJ504394)

STJ504394

GENERAL INFORMATION

Product Type	Proteins
Short Description	Myocilin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4653
Gene Symbol	MYOC
Uniprot ID	MYOC_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on Myocilin protein.
Immunogen Region	
Specificity	
Immunogen Sequence	