

## Myocilin protein (STJ504395)

STJ504395

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Myocilin-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	WB: IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	4653
<b>Gene Symbol</b>	MYOC
<b>Uniprot ID</b>	MYOC_HUMAN
<b>Immunogen</b>	Synthetic peptide corresponding to unique amino acid sequence on Myocilin protein.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	