

MyoD1 protein (STJ504397)

STJ504397

GENERAL INFORMATION

Product Type	Proteins
Short Description	MyoD1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4654
Gene Symbol	MYOD1
Uniprot ID	MYOD1_HUMAN
Immunogen	Synthetic peptide taken within amino acid sequence 280-330 on human MyoD1 protein.
Immunogen Region	280-330
Specificity	
Immunogen Sequence	