

Nicotinic Acetylcholine Receptor alpha 7 protein (STJ504399)

STJ504399

GENERAL INFORMATION

Product Type	Proteins
Short Description	Nicotinic Acetylcholine Receptor alpha 7-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25302
Gene Symbol	Chrna7
Uniprot ID	ACHA7_RAT
Immunogen	Synthetic peptide taken within amino acid region 25-75 on mouse GTPCHI protein.
Immunogen Region	25-75
Specificity	
Immunogen Sequence	