

PCSK9 protein (STJ504402)

STJ504402

GENERAL INFORMATION

Product Type	Proteins
Short Description	PCSK9-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100102
Gene Symbol	Pcsk9
Uniprot ID	PCSK9_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 30-60 on the NARC1 protein.
Immunogen Region	30-60
Specificity	
Immunogen Sequence	