

NME2 + NM23A protein (STJ504404)

STJ504404

GENERAL INFORMATION

Product Type	Proteins
Short Description	NME2 + NM23A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4830
Gene Symbol	NME1
Uniprot ID	NDKA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-150 on NM23A protein selective for both NDKA and NDKB.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	