

KIRREL protein (STJ504407)

STJ504407

GENERAL INFORMATION

Product Type	Proteins
Short Description	KIRREL-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	170643
Gene Symbol	Kirrel1
Uniprot ID	KIRRR1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 650-700 on mouse Neph1 protein. Binds to extracellular portion.
Immunogen Region	650-700
Specificity	
Immunogen Sequence	