

## CD10 protein (STJ504413)

STJ504413

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	CD10-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	4311
<b>Gene Symbol</b>	MME
<b>Uniprot ID</b>	NEP_HUMAN
<b>Immunogen</b>	Multi-epitope synthetic peptides were used: a mixture of synthetic peptides used for NEPR-101AP, NEPR-112AP, NEPR-121AP, NEPR-131AP, NEPR-141AP.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	