

NKG2D protein (STJ504417)

STJ504417

GENERAL INFORMATION

Product Type	Proteins
Short Description	NKG2D-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24934
Gene Symbol	Klrk1
Uniprot ID	NKG2D_RAT
Immunogen	Synthetic peptide taken within amino acid region 180-230 on rat NKG2D protein.
Immunogen Region	180-230
Specificity	
Immunogen Sequence	