

nmt55/p54nrb protein (STJ504421)
STJ504421**GENERAL INFORMATION**

Product Type	Proteins
Short Description	nmt55/p54nrb-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4841
Gene Symbol	NONO
Uniprot ID	NONO_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 430-471 on human NonO protein.
Immunogen Region	430-471
Specificity	
Immunogen Sequence	