

Niemann Pick C1 protein (STJ504424)

STJ504424

GENERAL INFORMATION

Product Type	Proteins
Short Description	Niemann Pick C1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:100 ICC: 1:100 IF: 1:100
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4864
Gene Symbol	NPC1
Uniprot ID	NPC1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1228-1278 on human Niemann-Pick C1 protein.
Immunogen Region	1228-1278
Specificity	
Immunogen Sequence	