

Nurr1 protein (STJ504433)

STJ504433

GENERAL INFORMATION

Product Type	Proteins
Short Description	Nurr1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4929
Gene Symbol	NR4A2
Uniprot ID	NR4A2_HUMAN
Immunogen	Synthetic peptide corresponding to amino acid region 566-591. The selected peptide was post-synthetically modified to achieve highest antigenicity before coupling to carrier protein using hetero bifunctional cross linker. Sequence: CTQGLQRIFYLKLEDLVP
Immunogen Region	566-591
Specificity	
Immunogen Sequence	