

Olig2 protein (STJ504436)

STJ504436

GENERAL INFORMATION

Product Type	Proteins
Short Description	Olig2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10215
Gene Symbol	OLIG2
Uniprot ID	OLIG2_HUMAN
Immunogen	Synthetic peptide corresponding to Human Olig2 amino acids 136-156. The selected peptide was post-synthetically modified to achieve highest antigenicity before coupling to a carrier protein using hetero-bifunctional cross linker for immunogen prepara
Immunogen Region	136-156
Specificity	
Immunogen Sequence	